

# EPU MINI Technical Specification Sheet



DECTRIS reserves the right to replace components with equivalent or better parts based on the availability at our suppliers.

Feature	Details
Processor	Intel Xeon E5-2630 v3 Processor <ul style="list-style-type: none"> <li>Processor Frequency: 2.4 GHz</li> <li>Cores/Threads: 8C/16T</li> </ul>
Usable Buffer Size	340 GB <sup>1</sup>
Network Interface	2 x 10 GbE DA/SFP+ <sup>2</sup>
Network Cables	SFP+ to SFP+, 10 GbE, Copper Twinax, 3 m (DCU-EPU direct attach copper cable)
File Transfer Rate	> 700 MB/s <sup>3</sup>

<sup>1</sup> Ramdisk (tmpfs)

<sup>2</sup> 1 x 10 GbE DA/SFP+ available, 1 x 10 GbE used for connection to the DCU

<sup>3</sup> The recommended size per file amounts to 1 GB in order to achieve the specified transfer rate.

Created	Reviewed	Released
SK	AM	SB